



BelCare™

Every Beltone hearing aid comes with a partner included who will be right beside you every step of the way on your journey to better hearing. You also receive BelCare, our comprehensive lifetime service and care program that includes free annual hearing evaluations, free hearing aid inspections and cleanings for as long as you own your hearing aids, patient care phoneline, one year limited manufacturer's warranty and more.

Please visit [Beltone.com](https://www.Beltone.com) for more information.



Beltone has won many awards including the 2020 Big Innovation, 2020 Edison Award and Fast Company World Changing Ideas 2020 for breakthrough technology allowing direct streaming to Android operating system devices and preserving battery life via Bluetooth low energy, and the 2020 Silver Stevie Award for Front-Line Customer Service Team of the Year.

© 2020 GN Hearing Care Corporation. All rights reserved. Beltone is a trademark of GN Hearing Care Corporation. Apple, the Apple logo, iPhone, iPad, iPod touch and Apple Watch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Google, Android, Google Pixel, Google Play and the Google Play logo are trademarks of Google Inc. Dolby and the double-D symbol are registered trademarks of Dolby Laboratories. The Bluetooth word mark and logos are registered trademarks owned by Bluetooth SIG, Inc.

Beltone Canada

2 East Beaver Creek Road, Building 3
Richmond Hill, ON L4B 2N3
Canada
1-800-387-3744
[beltone.ca](https://www.beltone.ca)



CABT101-20.08-Rev.A




IMAGINE

a hearing solution tailored to you

Beltone **Imagine™**





The all-new Beltone Imagine

Experience a revolutionary hearing aid that's raising the standard for hearing solutions. After years of research and development, we created a device that's uniquely tailored to you. Our groundbreaking hearing aid technology uses your individual ear shape to deliver sound as nature intended. Now you'll feel empowered and confident in any environment knowing you'll have premium sound quality, a full understanding of your surroundings, and the best speech recognition – even in the most challenging environments.

Along with the exceptional sound comes Beltone's tailored care: our promise to be right beside you every step of the way on your journey to better hearing.

Introducing
the world's first
truly tailored
hearing experience



made possible by **M&RIE***.
(Microphone & Receiver-In-Ear)

90%
of listeners prefer M&RIE for
natural sound quality**

10%
M&RIE increases the ability to
locate sounds by

Studies prove
M&RIE delivers sound
exactly as a
unique ear would***

*Full-featured hearing aid with M&RIE, microphone and receiver in ear, plus two directional microphones. Groth J. (2020). An innovative RIE with microphone in the ear lets users "hear with their own ears". Beltone white paper.

** Compared to pinna compensation and omnidirectionality using traditional microphones.

*** In listening test with normal hearing participants.



Tailored
Sound



Tailored
Connectivity



Tailored
Care



Tailored
Design

Hear as **nature** **intended**

Introducing Beltone Imagine, the world's first Receiver-In-Ear (RIE) hearing aid that uses the shape and contours of your unique ear to provide sound tailored exactly to you.

We challenged ourselves to take our hearing innovations to the next level and created something truly groundbreaking with Beltone Imagine. Our new revolutionary hearing technology, that we call M&RIE (Microphone & Receiver-in-Ear), filters and delivers the most natural sound yet. M&RIE is an RIE hearing aid unlike any other. It works with your unique ear shape to collect sound just as nature intended.

For the very first time with an RIE hearing aid, we're now able to put a microphone directly in your ear, instead of behind it, using your unique ear contours the way nature intended. This brings you premium sound quality, a full understanding of your surroundings, and even better speech recognition than before – ensuring you'll always feel confident about your hearing, even in difficult listening environments.

When you wear Beltone Imagine, you can focus on what you want to hear, without losing out on the other sounds around you. We created all of this so it's even easier for you to stay engaged and enjoy conversations in social settings.



7



6

Make your day easy and enjoyable

Made to work in any situation

You'll feel more involved and comfortable in your social settings - and here's why: you're always in control of your hearing with Beltone Imagine. Thanks to our innovative CrossLink Directionality 3 technology, you'll enjoy having a complete understanding of your surroundings and be able to focus on the sounds that matter the most. The technology continuously monitors every sound in your environment, so you can quickly and easily choose which sound to listen to. CrossLink Directionality 3 also provides excellent speech understanding, without any disorientation, because it automatically adjusts to your situation.



Zero-in on the conversation

No matter how complex the environment, you'll always have the best possible speech clarity with our brand-new Ultra Focus feature. When you're in demanding situations, like a bustling restaurant or an outdoor event, and you want to fully zero-in on just the sounds in front of you, simply turn on the Ultra Focus program – via your hearing aid or the Beltone HearMax™ app on your smartphone. When Ultra Focus mode is on, both the right and left hearing aids are focused on the sound directly ahead of you, eliminating peripheral noise.

Now you can enjoy the ease of always having the best conversations in any location.

Made for
iPhone | iPad | iPod

Works with
android

Stay engaged with your favorites via our award-winning direct audio streaming from iOS and Android™ devices

Beltone Imagine provides the best streaming capabilities with Bluetooth® Low Energy for iOS devices and a growing range of Android devices. With this premium technology, you'll receive crystal-clear audio directly into your hearing aids.



You'll feel empowered and free with Beltone Imagine Bluetooth Low-Energy compatibility and its ability to stream all your favorite services directly to your hearing aids. Enjoy streaming phone calls, directions, internet radio, music, podcasts, and audiobooks. Our upgraded technology also includes mix-in audio streaming that eliminates connection delays. Now there is no need to worry about missing the first words of an important phone call or not hearing key driving directions.



Beltone has won many awards including the 2020 Big Innovation, 2020 Edison Award and Fast Company World Changing Ideas 2020 for breakthrough technology allowing direct streaming to Android operating system devices and preserving battery life via Bluetooth low energy.

Countless ways to connect

With the Beltone HearMax app you can easily control volume and adjust programs, noise reduction, speech focus and more, directly from your iPhone, iPad, iPod touch, Apple Watch, or Android devices. There's also a Find My Hearing Aid feature built right into the app, just in case. Simply download the app for free from the App Store or Google Play and start managing any sound setting at your own convenience.



Stay connected

The multiple connectivity features of Beltone Imagine keep you connected to all your favorite entertainment and communication technologies. Choose from a range of wireless accessories that offer extra ways to enjoy a personalized listening experience.



Phone Link 2

Remote Control 2

TV Link 2

Remote Control

myPAL Pro

myPAL Micro

Beltone TV Link 2

Turn up the volume by streaming sound from your TV directly to your Beltone Imagine hearing aids.

Beltone Remote Control and Remote Control 2

With one of these controls in your hand, you can adjust volume and change hearing aid programs for TV, radio, computer, etc. quickly and easily.

Beltone Phone Link 2

Enjoy hands-free phone calls and stream conversations directly to your hearing aids via Bluetooth services for easy, seamless phone calls.

Beltone myPAL Pro and myPAL Micro

Attach the devices to friends, colleagues, family, etc. for the clearest sound and speech understanding, no matter the location.


Care tailored to your needs



Whether you're beginning your hearing journey or upgrading your hearing aids, we're here for you – every step of the way. It's important to us that you get the most out of your hearing aids, so you can get the most out of life! Be it in-store meetings or from the comfort of your home via Beltone Remote Care Live, we'll be right beside you, ensuring your hearing aids always work exactly how you imagined.

Enjoy hearing care from anywhere


Get care from the comfort of your home with Beltone Remote Care Live. Our latest video chat technology lets you and your hearing care professional easily meet face-to-face, via the Beltone HearMax app, to quickly resolve any issues and remotely fine-tune your hearing aids. Now whenever you're in need, you can get customized care – even from afar.



How Beltone Remote Care Live works




1. Your hearing care professional (HCP) initiates a video call from their office at the time you both agreed upon.



2. You simply accept the call.



3. Your HCP makes real-time fine-tuning adjustments to your hearing aids according to your needs.



4. Your hearing aids are fine-tuned to ensure you won't miss out on what's most important to you!

A man with a beard and a woman with short grey hair are smiling at each other in a craft store. The man is holding a bouquet of white and pink flowers. The woman is holding a bouquet of orange flowers. In the background, there are shelves with various colorful ribbons and potted plants.

A design
taylor-made
for you

The Beltone Imagine hearing solution was created with your needs in mind. The modern styling of our brand-new design has just one touchpoint on the back of the ear sitting lower for greater discretion. Aids sit comfortably in your ears, making them easy to handle and even easier to wear. With three different models, two rechargeability choices and a wide variety of color choices to blend with hair and skin tones, you can customize to your preferences and wear your Beltone Imagine hearing aids with confidence.

Style that makes a difference



Designed to suit you perfectly

Beltone Imagine was made with you in mind. The new design was based on years of feedback from our patients that they want their hearing aids to be comfortable, discreet and reliable so they can easily move about their day. The lightweight, ergonomic design is available in three models with a choice of rechargeable or non-rechargeable batteries. You choose what best fits your lifestyle.

Protection and durability

The all new Beltone Imagine is even more durable than previous hearing aid models. Every Beltone Imagine comes with an IP68 classification which means it is a completely enclosed model that can withstand dust and water.

Choose your favorite color

Because no two people are the same, having a diverse color palette is important. Beltone Imagine comes in an array of colors, so you can choose the color that fits your hair, skin and style best. Select from Sparkling Silver, Champagne, Gold, Bronze, Warm Grey, Graphite, Espresso, and Deep Black.



Stay charged so you never miss out

If you've chosen rechargeable Beltone Imagine models, recharging is simple - place your hearing aids in the charger overnight and they are ready again in the morning. The is extremely durable – even after 4 years of use it will retain 80% of its capacity.

Standard Charger

Your Beltone Imagine hearing aids (63DRWC models only) come with a sleek portable charger that also serves as a carrying case. You will find using the charger to be amazingly simple.



- Charges device when connected to a power outlet
- Charges through contactless charging
- Back LED light displays the charger is connected to a power source

Premium Charger

An optional premium charger is also available with Beltone Imagine hearing aids (63DRWC models only). The charger itself has an internal battery with an impressive capacity to give you 3 full charges.



- Provides up to three full charges due to its onboard battery
- Charges through contactless charging
- Comes with two LED lights:
 - Back light displays status of the charger
 - Front light displays status of the hearing aid

Battery life

30 hours

16 hours

8 hours

Charging required

3 hours

1 hour

30 minutes