

Inspiring a **new** generation

Beltone **Serene**™

Let your clients share their wisdom
with tiny rechargeable hearing aids
for noisy everyday life



A new generation of hearing

Wisdom is shared and heard through life's stories and experiences, bringing generations together.

Help more clients share their wisdom and be inspired with Beltone Serene™ – tiny hearing aids for noisy everyday life.

Powered by Organic Hearing

Our unique approach to creating hearing solutions that sound natural, feel natural and connect your clients naturally to the world.

Build satisfaction from day 1

Just as wisdom is shared, so are knowledge and know-how. Beltone Serene is based on all the experience and expertise we have, combined into our tiniest-yet rechargeable hearing aids. Our data shows that the sound processing technology in Beltone Serene results in higher acceptance, increased user satisfaction, and top ratings by users for sound quality and speech-in-noise performance.¹


With Beltone Serene, there's never been a better way to give the best possible client satisfaction.



1. Compared to legacy technology. End-user NPS rating, data on file

Tiny yet powerful

Having the microRIE in your portfolio can help increase conversion. It's the best of both worlds. The tiny rechargeable design with all-day power is filled with our latest sound technology for hearing in noise, so your clients don't need to compromise on features or experience.



Fitting options for all needs

New program for first-time users

The new experience profile 'First Time User Onboarding' is ideal for initial sound comfort.

New 'Hear in Noise' program

Give your clients a great hearing experience in noise with a program tailored for just the times when they need extra help.

New solution for single-sided deafness

Offer a superb wireless CROS/biCROS solution with greater audibility all around.

Improved connectivity

Prepared for Bluetooth® Auracast™

Enabling users to access broadcast audio directly in public settings – such as Airports and Stadiums, and to any compatible device.

Easy calls with a double-tap on the ear

TapControl lets your clients answer calls easily with a double-tap on the ear or one of the hearing aids.

Beltone HearMax app

Use the Beltone HearMax™ app to empower your clients to take control of their hearing – connect to TV-Streamer+ and access the new 'Hear in Noise' program.

New TV-Streamer+

The new TV-Streamer+ gives your clients immediate speech clarity and audio at their preferred volume.



Your style options

IP68-rated all-weatherproof designs, created for mild to profound hearing loss and available in stylish colours that complement hair and skin tones.



Key features overview

Features	SER17	SER9	SER6	SER4
Sound Quality				
 WARP bands	★★★★ 17	★★★★ 14	★★☆☆ 12	★★☆☆ 12
 Environment recognition	★★★★ (Smart Gain Pro II)	★★☆☆ (Smart Gain Pro)	★★☆☆ (Smart Gain Pro)	n/a
 Directionality	★★★★ (CrossLink Directionality 4)	★★☆☆ (CrossLink Directionality 3)	★★☆☆ (CrossLink Directionality 2)	★★☆☆ (CrossLink Directionality)
 Speech intelligibility in noise	★★★★ (Ultra Focus 2)	★★☆☆ (Ultra Focus)	★★☆☆ (Adaptive Directionality)	★★☆☆ (Adaptive Directionality)
 Noise Reduction	★★★★ (5 settings)	★★☆☆ (3 settings)	★★☆☆ (2 settings)	★★☆☆ (On / Off)
 Comfort	★★★★ (3 settings)	★★★★ (3 settings)	★★☆☆ (On / Off)	★★☆☆ (On / Off)
 Feedback Management	★★★★ (Feedback Eraser 2)	★★★★ (Feedback Eraser 2)	★★★★ (Feedback Eraser 2)	★★★★ (Feedback Eraser 2)
 Low frequency boost	★★★★ (3 settings)	★★☆☆ (2 settings)	★★☆☆ (On / Off)	★★☆☆ (On / Off)
Connectivity				
 Direct audio streaming	★★★★ For iPhone, iPad, iPod touch and Android™ devices. Hands-free calls for iPhone and iPad*			
 Public connectivity	★★★★ Future ready with Auracast for new broadcast opportunities,telecoil built into RIE 13 battery style and in MyPAL Pro			
 Accessories	★★★★ TV-Streamer+, TV Link 2, Phone Link 2, MyPAL Pro, MyPAL Micro, Remote Control, Beltone HearMax App			
Fitting Flexibility				
 Wireless CROS/BiCROS	★★★★			
 Gain handles	17	14	12	12
 Programs	4	4	4	4
 Fitting options	yes	yes	yes	yes

*Hands-free calls are compatible with iPhone 11 or later, iPad Pro 12.9-inch (5th generation), iPad Pro 11-inch (3rd generation), iPad Air (4th generation), iPad mini (6th generation) or later, with software updates iOS 15.3 and iPadOS 15.3 or later.

Worldwide Headquarters

Beltone A/S
Lautrupbjerg 7
DK-2750 Ballerup
Denmark
Tel.: + 45 45 75 11 11
beltone.com

Australia

GN Hearing Australia Pty Ltd
Gate C, 19-25 Khartoum Rd
Macquarie Technology Park
Macquarie Park NSW 2113
Tel.: (free) 1800 658 955
beltone.com/en-au

New Zealand

GN Hearing New Zealand
Limited
Ground Floor, North Entrance
4 Fred Thomas Drive
Takapuna, Auckland, 0622
Tel.: (free) 0800 900 126
beltone.com/en-au

Contact us
or visit
[beltone.com](https://www.beltone.com)
to learn more



© 2023 GN Hearing A/S. All rights reserved. Beltone is a trademark of GN Hearing A/S. Apple, the Apple logo, iPhone, iPad and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Android, Google Play and the Google Play logo are trademarks of Google LLC. The Bluetooth word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The Auracast word mark and logos are trademarks owned by Bluetooth SIG, Inc.