



THERE'S NEVER BEEN  
A BETTER MOMENT

WIDEX **MOMENT**<sup>™</sup>

**WIDEX**

SOUND LIKE NO OTHER



## WHAT MAKES A GREAT HEARING AID?

Designed to help you participate more fully

Promotes hearing well to keep your mind sharp

Produces natural sound to make hearing easier

Provides flexible hearing in any environment

Customized solution to set you up for success



# HEARING LIKE YOU REMEMBER

Until now, even the best hearing aids sounded like hearing aids. That's all changed with the enhanced Widex Moment with PureSound™ technology. Unlike other devices, Moment processes sound with lightening speed, so there is none of the delay that causes that annoying “tinny” artificial sound. Just pure, natural sound so you hear every moment, with sound that over 90% of wearers are satisfied with.<sup>1</sup>



“WHEN I AM WEARING  
MY MOMENT HEARING AIDS,  
I FORGET I HAVE THEM IN.”  
—BARB G

<sup>1</sup> Laura Winther Balling, Oliver Townend, Dana Helmink (2021). Sound Quality for All – The Benefits of Ultra-Fast Signal Processing. To appear in Hearing Review.

# WIDEX SOUNDRELAX SOOTHING SOUNDS AIMS TO HELP YOU REDUCE STRESS

In today's busy world, we are all becoming increasingly aware of the effects stress can have on our minds and bodies. This is why finding new and innovative ways to overcome anxiety and to relax has never been more important.

**Did you know our Widex Moment hearing aids can help you regain your inner peace?**

Widex SoundRelax™ provides a library of sounds to choose from to help relaxation, well-being, and concentration. This unique feature is designed to help all hearing aid wears — with or without tinnitus — deal with anxiety, soothe the mind and boost concentration. SoundRelax, when activated in your hearing aid, is accessible in the Widex Moment app.



Download the Widex Moment app  
to get started

\* The Widex MOMENT app works with all 2.4 GHz compatible Widex MOMENT hearing aid models. Use Widex Tonelink app for all other models.





## MORE WAYS TO CONNECT

Widex Moment offers direct streaming from iPhone® and compatible Android™ devices so you can stream high-quality calls, video, and music directly to your hearing aids.

With the Moment app, you have full control of your hearing aids, so it's never been easier to fine-tune your hearing to get exactly the sound you prefer — in the moments that are important to you.

### **BETTER HEARING MADE EASY**

With MyGuide, a feature in the Moment app designed to help you get the most out of your hearing aids, find step-by-step instructions, videos, and short courses that will help you with everything from turning your hearing aids on to personalizing their sound.

### **SOUND MADE FOR YOU**

The MySound™ 2.0 section of the app uses the power of AI to help adjust your sound. You can choose to have the AI guide you through a series of sound comparisons to tailor the sound to your particular situation. Or you can select a sound profile that has been automatically generated from the preferences of other Moment users in similar situations all over the world!



Download the Widex Moment app  
to get started

## BUILT TO LAST

Widex Moment is designed for maximum moisture resistance and long life — because less corrosion means more durability.

The advanced, water-resistant nano coating and the fact that our rechargeable options have no battery door means there is less opportunity for moisture to sneak in — providing you with a listening experience free from disruptions.

## DISCREET RECHARGEABILITY

Convenience is discreet with Widex Moment. Although the sRIC R D model is the market's smallest rechargeable receiver-in-canal hearing aid, you can rest easy knowing you'll have all the power you need to hear every moment.

Take better care of your Widex Moment sRIC R D hearing aids with Widex Charge n Clean. The charger offers:

- four-hour charge for a full day use
- 30-minute charge for four hours of use
- Sanitizing function of earware
- Drying function

Widex Charge n Clean





## WIDEX MOMENT MODELS AND STYLES

Widex Moment is available in a variety of colors to match your skin tone, hair color, or personality. Below are a few of our most popular options. Your Hearing Care Professional will help you find the perfect fit.

### HEARING AID MODELS



sRIC R D

#### RECEIVER-IN-CANAL (RIC)

The case sits behind the ear with the receiver in the ear-tip instead of the case. This allows for a smaller design without sacrificing quality.

- sRIC R D is the smallest rechargeable hearing aid on the market with a modern design
- For mild to severe hearing loss



BTE R D

#### BEHIND-THE-EAR (BTE)

The case sits behind the ear with a hook over the top of the ear. Hearing aids behind the ear are the perfect choice if you want extra power in a stylish design.

- Rechargeable options available
- T-coil comes standard
- For mild to severe-to-profound hearing loss



CIC

#### IN-THE-EAR (ITE)

The custom hearing aid solution sits discreetly in your ear. This solution is available in a variety of models from the easy handling ITE to the ultra small completely-in-the canal hearing aid.

- Customized to fit your ear canal
- For mild to moderately severe hearing loss

# WIRELESS ACCESSORIES

Accessories help you effortlessly connect to audio equipment, remote controls, television and mobile or landline phones. Ask your Hearing Care Professional for more details.

## BLUETOOTH CONNECTIVITY



reddot award 2018  
winner

### SOUND ASSIST\*

Multi-functional device featuring remote mic, hands-free, t-coil and remote control

### TV PLAY™\*\*

Stylish and stable TV streaming

## ADDITIONAL WIRELESS CONNECTIVITY\*\*



### COM-DEX

For hands-free streaming



### COM-DEX REMOTE MIC

For improved audibility in noise



### RC-DEX

Your remote control



### PHONE-DEX 2\*\*\*

For land line phone streaming



### TV-DEX

For enjoying TV



### FM+DEX

For flexible streaming

\* For compatible styles sRIC R D, RIC 312 D, BTE R D & BTE 13 D

\*\* For all styles with WidexLink™ technology (excludes CIC-M).

\*\*\* Additional handset sold separately. Must have a PHONE-DEX 2 to use additional handset.







SOUND LIKE NO OTHER

Widex is a registered trademark of WSAUD A/S.  
CBR02M/0822

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play and Google Play logo are trademarks of Google Inc. Subject to change with no further notice.



[WWW.WIDEX.COM](http://WWW.WIDEX.COM)